

## SIDE BY SIDE COMPARISON

### Product Category: 0.9% Sodium Chloride

**Current Code:**  
**2B1322Q**

SOD CHLOR 0.9% INJ USP  
VIAFLEX PLASTIC 250ML

Barcode on Bag ✓

Graduation Marks ✓

**Administration Port Closure Type:**  
Pull off port protector (blue color)




**Import Code:**  
**A6C1322**

SOD CHLOR 0.9% INJ USP  
IVINA PLASTIC 250ML

Barcode on Bag ✗

Graduation Marks ✓

**Administration Port Closure Type:**  
Twist off port protector (white color)



Translated  
product label

Key Product Differences	Current Code	Import Code
Product Code Number	2B1322Q	A6C1322
IV Bag Material	PVC	nPVC/nDEHP
Volume Location	Below product description	Below product description
Lot and manufacturing location	Top left - Lot number, Top right - Expiration Date	Bottom left; listed in order: Lot Number, Manufacturing Date, Expiration Date
Medication Port Location	Left side - Standard	Right side – Standard
Needle Type for Puncture of Medication Port	19-22G	18G
Compatible Spike Port	Designed to conform with ISO8536-4	Designed to conform with ISO8536-4
Vial-Mate Compatibility	Yes	No

For up-to-date information on these products, please visit the Baxter Medical Education Site using the QR code below.

