

SIDE BY SIDE COMPARISON

Product Category:
Balanced Solutions

Current Code:
2B2544X

Barcode on Bag ✓

Graduation Marks ✓

Administration Port Closure Type:
Pull off port protector (blue color)



Import Code:
GCCE0324

Barcode on Bag ✗

Graduation Marks ✓

Administration Port Closure Type:
Twist off port protector (white color)



Key Product Differences	Current Code & Description	Import Code & Description
Product Code Number	2B2544X: Plasma-Lyte A Injection pH 7.4 (Multiple Electrolytes Injection, Type 1, USP) in Viaflex Plastic Container, 1000 mL	GCCE0324: Plasma-Lyte 148 (pH 7.4) Solution for Infusion in Viaflo Plastic Container, 1000mL (no barcode on primary label)
IV Bag Material	PVC	nPVC/nDEHP
Volume Location	Below product description	Upper right
Lot and manufacturing location	Top left - Lot Number, Top right - Expiration Date	Bottom left - Lot Number, Bottom right - Expiration Date
Medication Port Location	Left side - Standard	Left side - Standard
Needle Type for Puncture of Medication Port	19-22G	19-22G
Compatible Spike Port	Designed to conform with ISO8536-4	Designed to conform with ISO8536-4

For up-to-date information on these products, please visit the Baxter Medical Education Site using the QR code below.

