

SIDE BY SIDE COMPARISON

Product Category:

0.9% Sodium Chloride Irrigation

Current Code:

2B7127

0.9% SODIUM CHLORIDE IRRIGATION, USP,
3000 ML UROMATIC PLASTIC CONTAINER

Barcode on Bag ✓

Graduation Marks ✓

Administration Port Closure Type:
Pull off port protector (blue color)



Import Code:

3KB7127

0.9% SODIUM CHLORIDE IRRIGATION, USP,
3000 ML UROMATIC PLASTIC CONTAINER

Barcode on Bag ✗

Graduation Marks ✓

Administration Port Closure Type:
Pull off port protector (blue color)



Key Product Differences	Current Code & Description	Import Code & Description
Product Code Number	2B7127: 0.9% SODIUM CHLORIDE IRRIGATION, USP, 3000 ML UROMATIC PLASTIC CONTAINER	3KB7127: 0.9% SODIUM CHLORIDE IRRIGATION, USP, 3000 ML UROMATIC PLASTIC CONTAINER
IV Bag Material	PVC	PVC
Volume Location	Below product description	Top right
Lot and manufacturing location	Top Left - Lot Number, Top right - Expiration Date	Bottom Left - Lot Number, Bottom Right - Expiration Date
Medication Port Location	Left side - Standard	Left side - Standard
Needle Type for Puncture of Medication Port	19-22G	21G
Compatible Spike Port	Designed to conform with ISO8536-4	Designed to conform with ISO8536-4

For up-to-date information on these products, please visit the Baxter Medical Education Site using the QR code below.

