

SIDE BY SIDE COMPARISON
Product Category:
0.9% Sodium Chloride

Current Code:
2B1323Q

SOD CHLOR 0.9% INJ USP
 VIAFLEX PLASTIC 500ML

Barcode on Bag ✓

Graduation Marks ✓

Administration Port Closure Type:
 Pull off port protector (blue color)




Import Code:
A6C1323US

SOD CHLOR 0.9% INJ USP
 IVINA PLASTIC 500ML

Barcode on Bag ✓

Graduation Marks ✓

Administration Port Closure Type:
 Twist off port protector (white color)



Key Product Differences	Current Code	Import Code
Product Code Number	2B1323Q	A6C1323US
IV Bag Material	PVC	nPVC/nDEHP
Volume Location	Below product description	Below product description
Lot and manufacturing location	Top left - Lot number, Top right - Expiration Date	Bottom left; listed in order: Lot Number, Manufacturing Date, Expiration Date
Medication Port Location	Left side - Standard	Right side – Standard
Needle Type for Puncture of Medication Port	19-22G	18G
Compatible Spike Port	Designed to conform with ISO8536-4	Designed to conform with ISO8536-4

For up-to-date information on these products, please visit the Baxter Medical Education Site using the QR code below.

