

**SIDE BY SIDE COMPARISON**

**Product Category:**

**5% Dextrose and 0.9% Sodium Chloride**

**Current Code:**  
**2B1064X**

5% DEXTROSE AND 0.9% SODIUM CHLORIDE INJECTION, USP, 1000 ML VIAFLEX PLASTIC CONTAINER

Barcode on Bag ✓

Graduation Marks ✓

**Administration Port Closure Type:**  
Pull off port protector (blue color)




**Import Code:**  
**A6C1064US**

5% GLUCOSE AND 0.9% SODIUM CHLORIDE INJECTION, USP, 1000 ML IVINA PLASTIC CONTAINER

Barcode on Bag ✓

Graduation Marks ✓

**Administration Port Closure Type:**  
Twist off port protector (white color)



Key Product Differences	Current Code	Import Code
Product Code Number	2B1064X	A6C1064US
IV Bag Material	PVC	nPVC/nDEHP
Volume Location	Below product description	Below product description
Lot and manufacturing location	Top left - Lot number, Top right - Expiration Date	Bottom left; listed in order: Lot Number, Manufacturing Date, Expiration Date
Medication Port Location	Left side - Standard	Right side – Standard
Needle Type for Puncture of Medication Port	19-22G	18G
Compatible Spike Port	Designed to conform with ISO8536-4	Designed to conform with ISO8536-4

For up-to-date information on these products, please visit the Baxter Medical Education Site using the QR code below.

